



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Lenz two-way plug valve with PTFE plug, valve NS 14.5, bore 1.5 mm, capillary approach tubes

**€93.63**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Two-way plug valves with PTFE plug

With threaded valve lock and two oblique bores.

Tightness meets the requirements of DIN 12540-1.

With capillary attachment tubes.

## Additional Information

---

No.	LZ-0369
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2368
EAN	4058072301118
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	1.5 mm
Material tap	Borosilicate glass
Material chicks	Polytetrafluoroethylene (PTFE)
Standard ground	44695
Type cock / valve	Tap 2-way
Type connection	Glass spigot

