



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz two-way plug valve with PTFE plug, valve NS 14.5, bore 2.5 mm, capillary approach tubes

€93.63
plus VAT &
Shipping

Product Images



Description

Two-way plug valves with PTFE plug

With threaded valve lock and two oblique bores.

Tightness meets the requirements of DIN 12540-1.

With capillary attachment tubes.

Additional Information

No.	LZ-0370
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2369
EAN	4058072301125
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	2.5 mm
Material tap	Borosilicate glass
Material chicks	Polytetrafluoroethylene (PTFE)
Standard ground	44695
Type cock / valve	Tap 2-way
Type connection	Glass spigot

