




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz vacuum plug valve, two-way, right-angled, valve NS 14.5

€50.21
plus VAT &
Shipping

Product Images



Description

Vacuum plug valves with right-angled lug tubes and hollow glass plug.

Additional Information

No.	LZ-0411
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2411
EAN	4058072301538
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	2.5 mm
Material tap	Borosilicate glass
Material chicks	Borosilicate glass
Standard ground	44695
Type cock / valve	Tap 2-way
Type connection	Glass spigot

