



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz vacuum plug valve, two-way, right angle, valve NS 29

€99.47
plus VAT &
Shipping

Product Images



Description

Vacuum plug valves with right-angled lug tubes and hollow glass plug.

Additional Information

No.	LZ-0415
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2415
EAN	4058072301576
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
sterile	No
Bore	10 mm
Material tap	Borosilicate glass
Material chicks	Borosilicate glass
Standard ground	29
Type cock / valve	Tap 2-way
Type connection	Glass spigot

