



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® tube connector (PP) with LL male, f. 6 mm I.-Ø, 10 pcs./pack

€26.10
plus VAT &
Shipping

Product Images



Description

High-quality hose Luer connectors with olive and Luer-Lock cone (male)

With precisely machined olives, to which also fixed hoses can be adapted.

The olives expand the hose up to 1.5 times the inner diameter.

PU: 10 pcs./pack

Additional Information

| | |
|--------------------------|-----------------------|
| No. | 2-1897 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072058746 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (violet) |
| Material | Polypropylene (PP) |
| Gender Luer lock | male |
| Type Hose Luer Connector | Olive/ Luer lock male |
| for DM hose inside 1 | 6.0 mm |

