



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® 2-in-1 plug-in spout without shut-off, acetal, male part, 2x1.7 mm Ø

**€17.10**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

With the 2-in-1 quick coupling, two separate flow paths can be easily connected or disconnected with 1 click. The freely rotatable quick coupling prevents kinks in the hoses.

## Additional Information

---

No.	2-0830
Manufacturer (Brand)	neoLab
EAN	4058072090272
Color	White
Material	Acetal
Form	straight
Length	28.4 mm
TBST MAX	82 °C
TBST MIN	-40 °C
Clutches properties	with hose olive
MAX operating pressure	8.3 bar abs.
Gender coupling part	male
for DM hose inside 1	44682
Type coupling part	Coupling part male. with olive

