



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® quick coupling, 1/4", male, shut-off,
hose olive 6.4 mm**

**€10.90
plus VAT &
Shipping**

Product Images



Abbildung ähnlich

Description

The neoLab quick couplings with a nominal width of 1/4" are ideally suited as a low-cost variant for all standard applications. The matching female coupling part has a plastic release button. This is especially advantageous for applications for which metal parts are not suitable.

- No visible metal parts
- Prevents accidental disconnection
- No leakage, shorter installation length, faster assembly

Additional Information

| | |
|------------------------|--------------------------------|
| No. | 2-0248 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072165086 |
| Color | White |
| Material | Acetal |
| Form | straight |
| Length | 53 mm |
| TBST MAX | 82 °C |
| TBST MIN | -40 °C |
| Clutches properties | with hose olive with shut-off |
| MAX operating pressure | 8.3 bar abs. |
| Gender coupling part | male |
| for DM hose inside 1 | 6 mm |
| Type coupling part | Coupling part male. with olive |

