



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Quick coupling, 1/4", male, acetal, thread 1/4" with Viton O-ring

€3.88
plus VAT &
Shipping

Product Images



Abbildung ähnlich

Description

The neoLab quick couplings with a nominal width of 1/4" are ideally suited as a low-cost variant for all standard applications. The matching female coupling part has a plastic release button. This is especially advantageous for applications for which metal parts are not suitable.

- No visible metal parts
- Prevents accidental disconnection
- No leakage, shorter installation length, faster assembly

Additional Information

No.	2-0258
Manufacturer (Brand)	neoLab
EAN	4058072189211
Color	White
Material	Acetal
Form	straight
Length	35.3 mm
TBST MAX	82 °C
TBST MIN	-40 °C
Clutches properties	with hose olive
MAX operating pressure	8.3 bar abs.
Gender coupling part	male
for DM hose inside 1	6 mm
Type coupling part	Coupling part male. with olive

