



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® quick coupling, 1/4", male, shut-off,  
elbow, hose olive 9.5 mm**

**€16.20  
plus VAT &  
Shipping**

## Product Images



Abbildung ähnlich

## Description

---

Proven in thousands of applications worldwide, neoLab's 1/4" nominal size quick disconnects offer the most comprehensive selection of sizes and configurations. Injection molded from thermoplastic acetal or chemically resistant polypropylene, the quick couplings are rugged and can be used with most mild chemicals.

- Integrated Connectors
- "Click" when connecting
- Connect and disconnect with one hand
- No leakage, shorter installation length, faster assembly

## Additional Information

---

No.	2-0260
Manufacturer (Brand)	neoLab
EAN	4058072155582
Color	White
Material	Acetal
Form	angled
TBST MAX	82 °C
TBST MIN	-40 °C
Clutches properties	with hose olive
MAX operating pressure	8.3 bar abs.
Connection for counterpart	1/4"
Gender coupling part	male
for DM hose inside 1	6.4 mm
Type coupling part	Coupling part male. with olive

