



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® quick coupling, 1/4", male, shut-off, hose olive 6.4 mm

**€18.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Proven in thousands of applications worldwide, neoLab's 1/4" nominal size quick disconnects offer the most comprehensive selection of sizes and configurations. Injection molded from thermoplastic acetal or chemically resistant polypropylene, the quick couplings are rugged and can be used with most mild chemicals.

- Integrated Connectors
- "Click" when connecting
- Connect and disconnect with one hand
- No leakage, shorter installation length, faster assembly

## Additional Information

---

No.	2-0870
Manufacturer (Brand)	neoLab
EAN	4058072048259
Color	White
Material	Acetal
Form	straight
TBST MAX	82 °C
TBST MIN	-40 °C
Clutches properties	with hose olive with shut-off
Figure	6
MAX operating pressure	8.3 bar abs.
Connection for counterpart	1/4"
Gender coupling part	male
for DM hose inside 1	6.4 mm
Type coupling part	Coupling part male. with olive

