



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® quick coupling, 1/4", male, hose olive 6.4 mm

€6.31
plus VAT &
Shipping

Product Images



Description

Proven in thousands of applications worldwide, neoLab's 1/4" nominal size quick disconnects offer the most comprehensive selection of sizes and configurations. Injection molded from thermoplastic acetal or chemically resistant polypropylene, the quick couplings are rugged and can be used with most mild chemicals.

- Integrated Connectors
- "Click" when connecting
- Connect and disconnect with one hand
- No leakage, shorter installation length, faster assembly

Additional Information

No.	2-1871
Manufacturer (Brand)	neoLab
EAN	4058072058487
Transport temperature	Room temperature
Color	beige
Material	Polypropylene (PP)
Form	straight
TBST MAX	82 °C
TBST MIN	-40 °C
Clutches properties	with hose olive
Figure	7
MAX operating pressure	8.3 bar abs.
Connection for counterpart	1/4"
Gender coupling part	male
for DM hose inside 1	6.4 mm
Type coupling part	Coupling part male. with olive

