



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Quick coupling, 1/4", male, elbow, hose olive 6.4 mm

€14.60
plus VAT &
Shipping

Product Images



Description

Proven in thousands of applications worldwide, neoLab's 1/4" nominal size quick disconnects offer the most comprehensive selection of sizes and configurations. Injection molded from thermoplastic acetal or chemically resistant polypropylene, the quick couplings are rugged and can be used with most mild chemicals.

- Integrated Connectors
- "Click" when connecting
- Connect and disconnect with one hand
- No leakage, shorter installation length, faster assembly

Additional Information

No.	2-1874
Manufacturer (Brand)	neoLab
EAN	4058072058517
Transport temperature	Room temperature
Color	beige
Material	Polypropylene (PP)
Form	angled
TBST MAX	82 °C
TBST MIN	-40 °C
Clutches properties	with hose olive
Figure	8
Connection for counterpart	1/4"
MAX operating pressure	8.3 bar abs.
Gender coupling part	male
Type coupling part	Coupling part male. with olive
for DM hose inside 1	6.4 mm

