



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Quick coupling, 1/8", male, shut-off,
elbow, hose olive 3.2 mm**

**€21.40
plus VAT &
Shipping**

Product Images



Description

The male neoLab quick couplings with a nominal diameter of 1/8" cover a wide range of standard applications. With protective hose collar for additional protection of the 1.6 mm hose olive

- No leakage, shorter installation length, faster assembly
- Stainless steel release button
- "Click" when connecting with matching female counterpart

Additional Information

No.	2-0850
Manufacturer (Brand)	neoLab
EAN	4058072157524
Color	ivory
Material	Polypropylene (PP)
Form	angled
Length	27.6 mm
TBST MAX	71 °C
TBST MIN	0 °C
Clutches properties	with hose olive
MAX operating pressure	8.3 bar abs.
Connection for counterpart	1/8"
Gender coupling part	male
for DM hose inside 1	3.2 mm
Type coupling part	Coupling part male. with olive

