



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Quick coupling, 1/8", male, elbow, hose olive 6.4 mm

€14.30
plus VAT &
Shipping

Product Images



Description

The male neoLab quick couplings with a nominal diameter of 1/8" cover a wide range of standard applications.

- No leakage, shorter installation length, faster assembly
- Stainless steel release button
- "Click" when connecting with the matching female counterpart

Additional Information

No.	2-1867
Manufacturer (Brand)	neoLab
EAN	4058072058449
Transport temperature	Room temperature
Color	beige
Material	Polypropylene (PP)
Form	straight
TBST MAX	71 °C
TBST MIN	0 °C
Clutches properties	with hose olive
MAX operating pressure	8.3 bar abs.
Connection for counterpart	1/8"
Gender coupling part	male
for DM hose inside 1	6.4 mm
Type coupling part	Coupling part male. with olive

