



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® quick coupling, 3/8", male, shut-off, hose olive 9.5 mm

€21.70
plus VAT &
Shipping

Product Images



Description

The compact and lightweight neoLab quick couplings made of chemically resistant polypropylene with a nominal diameter of 3/8" replace the heavy and bulky standard brass couplings with ball locks in a wide range of applications. Ergonomic design and large release buttons on the female counterpart result in a coupling that is easy to handle and operate. The particularly efficient valve design ensures high flow and exceptionally low dead volume.

- optionally with hose olive or with male thread
- with shut-off

Additional Information

No.	2-0812
Manufacturer (Brand)	neoLab
EAN	4058072048099
Color	gray
Material	Polypropylene (PP)
Form	straight
TBST MAX	110 °C
TBST MIN	0 °C
Clutches properties	with hose olive
Figure	6
MAX operating pressure	7.2 bar abs.
Connection for counterpart	3/8"
Gender coupling part	male
for DM hose inside 1	9.5 mm
Type coupling part	Coupling part male. with olive

