



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® quick coupling, 9/32", male, shut-off, hose olive 6.4 mm

**€23.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

The **neoLab 9/32" nominal size quick couplings** have a high performance valve that provides greater flow capacity than any other coupling of this size.

Chemically resistant polypropylene is ideal for use in extreme conditions. Some feature an optional bulkhead fitting for tight mounting to tank walls and drums.

- Greater flow despite compact size

## Additional Information

No.	2-0804
Manufacturer (Brand)	neoLab
EAN	4058072048013
Color	gray
Material	Polypropylene (PP)
Form	straight
TBST MAX	110 °C
TBST MIN	0 °C
Clutches properties	with hose olive
Figure	5
MAX operating pressure	7.2 bar abs.
Connection for counterpart	9/32"
Gender coupling part	male
for DM hose inside 1	6.5 mm
Type coupling part	Coupling part male. with olive

