



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Quick coupling, 1/8", female, shut-off, hose olive 6.4 mm

**€19.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The robust neoLab quick couplers with a nominal diameter of 1/8" and a stainless steel release button with "click" function provide an even more stable hold when connecting hoses. Available in thermoplastic acetal or chemically resistant polypropylene.

- Ensures a secure connection
- No leakage, shorter installation length, faster assembly

## Additional Information

---

No.	2-1863
Manufacturer (Brand)	neoLab
EAN	4058072058401
Transport temperature	Room temperature
Color	beige
Material	Polypropylene (PP)
Form	straight
TBST MAX	71 °C
TBST MIN	0 °C
Clutches properties	with hose olive with shut-off
Figure	4
MAX operating pressure	8.3 bar abs.
Connection for counterpart	1/8"
Gender coupling part	female
for DM hose inside 1	6.4 mm
Type coupling part	Coupling part female. with olive

