



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Plug-in nozzle with olive, brass, NW 2.7 for 4 mm Ø

Special Price

€0.60 was

€0.90

plus VAT &
Shipping

Product Images



Description

One-way shut-off quick-connect couplings.
Versatile connection options are guaranteed.
Large passage with small construction dimensions and low
Pressure drop during coupling.
Material: Brass

Additional Information

No.	2-6292
Manufacturer (Brand)	neoLab
EAN	4058072073572
Transport temperature	Room temperature
Color	Gold
Material	Brass
Form	straight
Length	24 mm
TBST MAX	100 °C
TBST MIN	-20 °C
Clutches properties	with hose olive
MAX operating pressure	35 bar abs.
Gender coupling part	male
Type coupling part	Coupling part male. with olive

