




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® screw coupling with PTFE seal, GL 25

€31.16
plus VAT &
Shipping

Product Images



Description

Screw coupling made of chemically resistant PBT and integrated silicone (VMQ) seal... The screw coupling can be used to flexibly connect two glass threads. Ideal for removing solvents directly from storage vessels without the need for decanting.

Thread: GL 25

Additional Information

No.	SD-2723
Manufacturer (Brand)	DURAN
old neoLab article no.	E-2922
EAN	4058072363314
Transport temperature	Room temperature
Color	red white
Material	Polybutylene terephthalate (PBT)
Form	straight
Length	48 mm
Connection for counterpart	GL 25
Gender coupling part	female
Thread [connection]	GL 25
TBST MAX	180 °C
TBST MIN	-45 °C
Type coupling part	Screw coupling with GL union couplings

