



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® transition piece (PP) straight, f. 2.4-3.4 mm/3.2-4.2 mm, 10 pcs./pack

€25.50
plus VAT &
Shipping

Product Images



Description

Reducer, straight, made of natural-colored PP

Different olive size on both sides. Slightly transparent, **autoclavable**. Outstanding quality with constant inner diameter, causes uniform flow rates and tight connection. The olives are machined on the outside of the sealing cone so that solid tubing can be attached. They expand the hose up to 1.5 times its inner \varnothing .

In a bag of 10 pieces.

Additional Information

| | |
|-----------------------|-----------------------|
| No. | 2-0955 |
| Manufacturer (Brand) | neoLab |
| EAN | 4058072048471 |
| autoclavable at 121°C | Yes |
| Color | transparent (milky) |
| Material | Polypropylene (PP) |
| Form | straight |
| Symmetry | for different hose DM |
| Type hose connector | Reducer (straight) |
| for DM hose inside 1 | 2.4 to 3.4 mm |
| for DM hose inside 2 | 3.2 to 4.2 mm |

