



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® transition piece (PP) straight, f. 3.2-4.2
mm/6.4-7.4 mm, 10 pcs./pack**

**€25.50
plus VAT &
Shipping**

Product Images



Description

Reducer, straight, made of natural-colored PP

Different olive size on both sides. Slightly transparent, **autoclavable**. Outstanding quality with constant inner diameter, causes uniform flow rates and tight connection. The olives are machined on the outside of the sealing cone so that solid tubing can be attached. They expand the hose up to 1.5 times its inner \varnothing .

In a bag of 10 pieces.

Additional Information

No.	2-0961
Manufacturer (Brand)	neoLab
EAN	4058072048532
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Form	straight
Symmetry	for different hose DM
Type hose connector	Reducer (straight)
for DM hose inside 1	3.2 to 4.2 mm
for DM hose inside 2	6.4 to 7.4 mm

