



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz core grinding NS 45/40, 40x2,3x145 mm, 10 pcs/pack

**€115.70
plus VAT &
Shipping**

Product Images

Description

Ground taper cores made of DURAN® pipe. According to DIN 12249 Form B and in addition to the standard.

Additional Information

No.	LZ-0027
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2026
EAN	4058072297695
Transport temperature	Room temperature
diameter	40 mm
Material	Borosilicate glass
Length	145 mm
Standard ground core	45/40
standardized according to DIN standard (DIN)	DIN 12249

