



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Thread adapter borosilicate 3.3, GL 14, sleeve NS 14.6

€15.50
plus VAT &
Shipping

Product Images



Description

Thread adapters made of borosilicate glass 3.3 in various sizes. These are required to attach threaded vials to ground joint apparatus such as rotary evaporators. Size-dependent screw couplings (E-2920 ff) are required for this purpose.
GL 14, sleeve NS 14.5

Additional Information

No.	E-2930
Manufacturer (Brand)	neoLab
EAN	4058072086657
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Reducers Features	with standard ground joint
Standard ground sleeve	14/23
Type Connection Output	Screw connection GL 14
Type Connection Input	Standard ground

