



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® T-connector, 1 x 3.2-4.2 mm, 2 x 2.4-3.4 mm, PP, 10 pcs/pack

**€25.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Reducer, T-shape, made of natural-colored PP

**With different olive size.** Slightly transparent, **autoclavable**. Outstanding quality with constant inner diameter, causes uniform flow rates and tight connection. The olives are worked on the outside of the sealing cone in such a way that solid tubing can also be plugged on. They expand the hose up to 1.5 times its inner  $\varnothing$ .

**In a bag of 10 pieces.**

## Additional Information

---

No.	2-0974
Manufacturer (Brand)	neoLab
EAN	4058072048662
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Form	T-shaped
Symmetry	for different hose DM
Type hose connector	Reducer (T-shape)
for DM hose inside 1	2.4 to 3.4 mm
for DM hose inside 2	2.4 to 3.4 mm
for DM hose inside 3	3.2 to 4.2 mm

