



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® T-connector, 1 x 2.4-3.4 mm, 2 x 3.2-4.2 mm, PP, 10 pcs/pack

**€25.50
plus VAT &
Shipping**

Product Images



Description

Reducer, T-shape, made of natural-colored PP

With different olive size. Slightly transparent, **autoclavable**. Outstanding quality with constant inner diameter, ensures uniform flow rates and tight connection. The olives are worked on the outside of the sealing cone in such a way that solid tubing can also be plugged on. They expand the hose up to 1.5 times its inner \varnothing .

In a bag of 10 pieces.

Additional Information

No.	2-0979
Manufacturer (Brand)	neoLab
EAN	4058072048716
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polypropylene (PP)
Form	T-shaped
Symmetry	for different hose DM
Type hose connector	Reducer (T-shape)
for DM hose inside 1	2.4 to 3.4 mm
for DM hose inside 2	3.2 to 4.2 mm
for DM hose inside 3	3.2 to 4.2 mm

