



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® WILO clamp 15 mm

€6.85
plus VAT &
Shipping

Product Images



Description

Hose clamps or pinchcocks, open at the side due to folding mechanism, very sturdy design made of galvanized steel. With easy-grip, large knurled screw, easy and accurate adjustment.

Temperature resistant: up to 150°C

Additional Information

No.	DT-0611
Manufacturer (Brand)	neoLab
old neoLab article no.	1-3015
EAN	4058072228873
Transport temperature	Room temperature
Color	silver
Wide	4 mm
Length	57 mm
suitable for	Hoses
Material clamp	Steel galvanized
MAX span	15 mm
Span MIN	0 mm

