




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® hose clamp stainless, 20-32 mm clamping width

€4.34 As low as
€1.45
plus VAT &
Shipping

Product Images



Description

Stainless hose clamps with worm thread for secure and permanent fastening of hoses.
Clean and hygienic, acid resistant.
Hexagonal screw with slot and cross recess.

Additional Information

No.	DT-0635
Manufacturer (Brand)	neoLab
old neoLab article no.	2-3073
EAN	4058072228972
Transport temperature	Room temperature
Color	silver
Material	Stainless steel
Wide	9 mm
suitable for	Hoses
Clamps properties	with special chemical resistance
MAX span	32 mm
Span MIN	20 mm

