



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolabLine neoClips hose clamps size BB, 9.2 - 10.3 mm, 10 pcs./pack

€12.10
plus VAT &
Shipping

Product Images



Description

neoClips are absolutely rustproof, unbreakable and **extremely heat-resistant**.

Even chemicals and other sharp substances cannot harm them.

With neoClips, hoses, tubes, etc. can be fastened **absolutely tightly and reliably**. The clips are easy to mount. By pressing the two ends together, the teeth interlock and snap into place. The clip can also be released at any time.

PU: 10 pcs./pack.

Technical data

pressure resistant up to 14 bar

UK listed and FDA certified

Additional Information

No.	1-0053
Manufacturer (Brand)	neoLabLine
EAN	4058072031022
Color	Black
Material	Polyamide 6.6 (PA 6.6. Nylon)
suitable for	Hoses
TBST MAX	125 °C
TBST MIN	-40 °C
Clamps properties	with special chemical resistance
Size hose clamp	BB
MAX span	10.3 mm
Span MIN	9.2 mm

