



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolabLine neoClips hose clamps size H, 18.1 – 20.6 mm, 10 pcs./pack

**€13.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

**neoClips** are absolutely rustproof, unbreakable and **extremely heat-resistant**.

Even chemicals and other sharp substances cannot harm them.

With neoClips, hoses, tubes, etc. can be fastened **absolutely tightly and reliably**. The clips are easy to install. By pressing the two ends together, the teeth interlock and snap into place. The clip can also be released at any time.

**PU:** 10 pcs./pack.

### Technical data

pressure resistant up to 14 bar

UK listed and FDA certified

## Additional Information

---

No.	1-0060
Manufacturer (Brand)	neoLabLine
EAN	4058072031091
Color	Black
Material	Polyamide 6.6 (PA 6.6. Nylon)
suitable for	Hoses
TBST MAX	125 °C
TBST MIN	-40 °C
Clamps properties	with special chemical resistance
Size hose clamp	H
MAX span	20.6 mm
Span MIN	18.1 mm

