



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Norprene tubing, 6.4 x 9.6 mm, 15
m/roll**

**€87.30
plus VAT &
Shipping**

Product Images



Description

Thermoplastic elastomer hose, black.

Suitable for peristaltic pumps. Very long service life. Resistant to many aggressive chemicals. Heat and ozone resistant, does not become brittle.

PU: 15 m/roll

Additional Information

No.	7-9249
Manufacturer (Brand)	neolab
EAN	4058072049874
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Black
Material	Elastomer (Norprene)
Form	smooth straight
Length	15 m
DM outside	9.6 mm
DM inside	6.4 mm
TBST MAX	135 °C
TBST MIN	-60 °C
Wall thickness hose	1.6 mm
Hoses properties	with special chemical resistance suitable for foodstuffs
Application hose	for peristaltic pumps
Hardness	61 shore A
Article meets criteria according to	FDA
for medium	Liquids Aqueous solutions

