



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Norprene tubing, 12.7 x 17.5 mm, 15 m/roll

€279.00
plus VAT &
Shipping

Product Images



Description

Thermoplastic elastomer hose, black.

Suitable for peristaltic pumps. Very long service life. Resistant to many aggressive chemicals. Heat and ozone resistant, does not become brittle.

PU: 15 m/roll

Additional Information

No.	7-9252
Manufacturer (Brand)	neolab
EAN	4058072049904
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	Black
Material	Elastomer (Norprene)
Form	smooth straight
Length	15 m
DM outside	17.5 mm
DM inside	12.7 mm
TBST MAX	135 °C
TBST MIN	-60 °C
Wall thickness hose	2.4 mm
Hoses properties	with special chemical resistance suitable for foodstuffs
Application hose	for peristaltic pumps
Hardness	61 shore A
Article meets criteria according to	FDA
for medium	Liquids Aqueous solutions

