



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® PharMed tubing for cell culture, 3.2 x 6.4 mm, 7.5 m/roll

€174.20
plus VAT &
Shipping

Product Images



Description

Hose made of thermoplastic elastomer.

Beige, non-toxic. Less gas permeable than silicone. Long service life in peristaltic pumps. Ozone, UV, chemical resistant.

PU: 7.5 m/roll

Additional Information

No.	7-9260
Manufacturer (Brand)	neoLab
EAN	4058072049911
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	beige
Material	Elastoplastic (TPE)
Form	smooth straight
Length	7.5 m
DM outside	6.4 mm
DM inside	3.2 mm
TBST MAX	135 °C
TBST MIN	-73 °C
Wall thickness hose	1.6 mm
Hoses properties	with special chemical resistance suitable for foodstuffs
Application hose	for peristaltic pumps
Hardness	61 shore A
Article meets criteria according to	FDA USP
for medium	Liquids Aqueous solutions

