



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® PharMed tubing for cell culture, 4.8 x 8.0 mm, 7.5 m/roll

€175.70  
plus VAT &  
Shipping

### Product Images



## Description

---

### Hose made of thermoplastic elastomer.

Beige, non-toxic. Less gas permeable than silicone. Long service life in peristaltic pumps. Ozone, UV, chemical resistant.

**PU:** 7.5 m/roll

## Additional Information

---

No.	7-9261
Manufacturer (Brand)	neoLab
EAN	4058072049928
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	beige
Material	Elastoplastic (TPE)
Form	smooth straight
Length	7.5 m
DM outside	8 mm
DM inside	4.8 mm
TBST MAX	135 °C
TBST MIN	-73 °C
Wall thickness hose	1.6 mm
Hoses properties	with special chemical resistance suitable for foodstuffs
Application hose	for peristaltic pumps
Hardness	61 shore A
Article meets criteria according to	FDA USP
for medium	Liquids Aqueous solutions

