



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Deutsch & Neumann silicone hose 1 x 1.5 mm, 1 m

**€2.35**  
plus VAT &  
Shipping

### Product Images



## Description

Standard silicone hose

Hardness: 55 Shore A

**Temperature resistant** at -60 to +180°C (short-term up to +260°C) also suitable for hose pumps.

Price/m

## Additional Information

No.	NM-0340
Manufacturer (Brand)	Deutsch & Neumann
old neoLab article no.	E-4193
EAN	4058072337445
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Silicone
Form	smooth straight
Length	1 m
DM outside	4 mm
DM inside	1 mm
TBST MAX	180 °C
TBST MIN	-60 °C
Wall thickness hose	1.5 mm
Hoses properties	food safe
Application hose	for general laboratory applications for peristaltic pumps
Hardness	55 shore A
Article meets criteria according to	BGA DAB FDA USP
for medium	Liquids Aqueous solutions Air/Gases/Vapors

