



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Deutsch & Neumann silicone hose 4 x 1.5 mm, 1 m

€65.75
plus VAT &
Shipping

Product Images



Description

Standard silicone hose

Hardness: 55 Shore A

Temperature resistant at -60 to +180°C (short-term up to +260°C) also suitable for hose pumps.

Price/m

Additional Information

No.	NM-0354
Manufacturer (Brand)	Deutsch & Neumann
old neoLab article no.	E-4197
EAN	4058072337483
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Silicone
sterile	No
Form	smooth straight
Length	1 m
DM outside	7 mm
DM inside	4 mm
TBST MAX	180 °C
TBST MIN	-60 °C
Wall thickness hose	1.5 mm
Hoses properties	food safe
Application hose	for general laboratory applications for peristaltic pumps
Hardness	55 shore A
Article meets criteria according to	BGA DAB FDA USP
for medium	Liquids Aqueous solutions Air/Gases/Vapors

