



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Silicone tubing, 6 x 10 mm, 5 m/roll

€39.00
plus VAT &
Shipping

Product Images



Description

Transparent to translucent tube.

Biologically neutral. Good physical properties, supple, elastic, flexible.

Density: 1.16.

Elongation at break: >300 %

Corresponds to Ph.Eur. (2nd ed.), DAB 10 and USP XXII Class VI (Biol. Tests).

Additional Information

No.	7-9272
Manufacturer (Brand)	neoLab
EAN	4058072049973
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Silicone
Form	smooth straight
Length	5 m
DM outside	10 mm
DM inside	6 mm
TBST MAX	200 °C
TBST MIN	-60 °C
Wall thickness hose	2 mm
Hoses properties	food safe
Application hose	for general laboratory applications for cell culture
Article meets criteria according to	FDA BGA EU guidelines
Hardness	60 shore A
for medium	Liquids Aqueous solutions Air/Gases/Vapors

