



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Socorex replacement aspiration tubing, uncut
per meter, for Calibrex 520/525/530 up to 10 ml**

**€26.78
plus VAT &
Shipping**

Product Images



Description

Standard tubing

Each ^{Calibrex™} dispenser is supplied with a suction tube.

Replacement tubing is available either cut to standard size or uncut, per meter.

7-7478: Metered/uncut

Additional Information

No.	SR-0202
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7478
EAN	4058072367206
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polytetrafluoroethylene (PTFE)
Length	1000 mm
DM inside	5 mm
Hoses properties	with special chemical resistance
Application hose	for Socorex Dispenser
for medium	Aqueous solutions

