



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex suction hose, uncut, per meter, for Calibrex 525/530 Dispensers 25/50/100 ml

€22.93
plus VAT &
Shipping

Product Images



Description

Standard tubing

Each ^{Calibrex™} dispenser is supplied with a suction tube.

Replacement tubing is available either cut to standard size or uncut, per meter.

7-7635: Metered/uncut

Additional Information

No.	SR-0254
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7635
EAN	4058072367725
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Perfluoroethylene propylene (FEP)
Length	1000 mm
Hoses properties	with special chemical resistance
Application hose	for Socorex Dispenser
DM inside	7 mm
for medium	Aqueous solutions

