



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Socorex Replacement Suction Tubing for Dosystm Self-Filling Syringes 0.3 to 10 mL, cut

**€7.90**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

### Medical grade silicone hoses

Suction tubing for syringes made of durable medical grade silicone.

Cut

Autoclavable at 121°C.

## Additional Information

---

No.	SR-0335
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7739
EAN	4058072368524
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Perfluoroethylene propylene (FEP)
Length	1000 mm
Hoses properties	with special chemical resistance
Application hose	for Socorex Dispenser
for medium	Aqueous solutions

