



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Telescopic Tubing, 25, 50 and 100 ml for Calibrex Dispensers 525/530, 195-345 mm

€19.66
plus VAT &
Shipping

Product Images



Description

Telescopic hoses
Simultaneous adaptation to any bottle size.

7-7637: Telescopic hose 195 to 345 mm

Additional Information

No.	SR-0255
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7637
EAN	4058072367732
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Perfluoroethylene propylene (FEP)
Hoses properties	with special chemical resistance
Application hose	for Socorex Dispenser
for medium	Aqueous solutions

