



neolab Migge GmbH
 Rischerstr. 7-9
 69123 Heidelberg
 Deutschland
 +49 (0)6221 /
 8442-44
<https://www.neolab.de>
 e

Umsatzsteuer-
 Identifikationsnummer
 :
 DE 143 450 657



neolab® Tygon peristaltic tubing SH 75, 6.4 x 9.6 mm, 15 m/roll

€179.00
 plus VAT &
 Shipping

Product Images



Description

A hose type of the extra class!

The exceptionally smooth inner wall prevents particle deposits.

Crystal clear, plasticizer-free, reduced risk of contamination.

Due to the low absorption capacity, the liquid remains unchanged.

REACH conformity.

Density: 0.9

PU: 15 m/roll

Additional Information

No.	2-4473
Manufacturer (Brand)	neoLab
EAN	4058072003562
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polyvinyl chloride (PVC. Tygon)
Form	smooth straight
Length	15 m
DM outside	9.6 mm
DM inside	6.4 mm
TBST MAX	54 °C
TBST MIN	-59 °C
Wall thickness hose	1.6 mm
Hoses properties	food safe
Application hose	for general laboratory applications for peristaltic pumps
Hardness	75 shore A
Article meets criteria according to	FDA
for medium	Liquids Aqueous solutions Air/Gases/Vapors

