



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Tygon peristaltic tubing SH 75, 9.5 x 12.7 mm, 15 m/roll

€206.00
plus VAT &
Shipping

Product Images



Description

A hose type of the extra class!

The exceptionally smooth inner wall prevents particle deposits.

Crystal clear, plasticizer-free, reduced risk of contamination.

Due to the low absorption capacity, the liquid remains unchanged.

REACH conformity.

Density: 0.9

PU: 15 m/roll

Additional Information

No.	2-4475
Manufacturer (Brand)	neoLab
EAN	4058072003586
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polyvinyl chloride (PVC, Tygon)
Form	smooth straight
Length	15 m
DM outside	12.7 mm
DM inside	9.5 mm
TBST MAX	54 °C
TBST MIN	-59 °C
Wall thickness hose	1.6 mm
Hoses properties	food safe
Application hose	for general laboratory applications for peristaltic pumps
Article meets criteria according to	FDA
Hardness	75 shore A
for medium	Liquids Aqueous solutions Air/Gases/Vapors

