




neolab Migge GmbH
 Rischerstr. 7-9
 69123 Heidelberg
 Deutschland
 +49 (0)6221 /
 8442-44
<https://www.neolab.de>
 e

Umsatzsteuer-
 Identifikationsnummer
 :
 DE 143 450 657



neolab® Tygon tubing S-hardness 55, 3.2 x 4.8 mm, 15 m/roll

€38.50
 plus VAT &
 Shipping

Product Images



Description

Ideal for general laboratory use.

High chemical resistance. Clear, soft and flexible. Lower oxygen permeability than rubber, aging resistant and non-oxidizing.

Quick and easy to plug onto fittings.

Additional Information

No.	2-4453
Manufacturer (Brand)	neoLab
EAN	4058072003449
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polyvinyl chloride (PVC. Tygon)
Form	smooth straight
Length	15 m
DM outside	4.8 mm
DM inside	3.2 mm
TBST MAX	74 °C
TBST MIN	-50 °C
Wall thickness hose	0.8 mm
Application hose	for general laboratory applications
Hardness	55 shore A
for medium	Liquids Aqueous solutions Air/Gases/Vapors

