



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Tygon tubing S-hardness 55, 4.8 x 6.4 mm, 15 m/roll

€46.70
plus VAT &
Shipping

Product Images



Description

Ideal for general laboratory use.

High chemical resistance. Clear, soft and flexible. Lower oxygen permeability than rubber, aging resistant and non-oxidizing.

Quick and easy to plug onto fittings.

Additional Information

No.	2-4456
Manufacturer (Brand)	neoLab
EAN	4058072003470
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polyvinyl chloride (PVC. Tygon)
Form	smooth straight
Length	15 m
DM outside	6.4 mm
DM inside	4.8 mm
TBST MAX	74 °C
TBST MIN	-50 °C
Wall thickness hose	0.8 mm
Application hose	for general laboratory applications
Hardness	55 shore A
for medium	Liquids Aqueous solutions Air/Gases/Vapors

