



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Tygon tubing S-hardness 55, 9.5 x 12.7 mm, 15 m/roll

**€182.70**  
plus VAT &  
Shipping

### Product Images



## Description

---

Ideal for general laboratory use.

High chemical resistance. Clear, soft and flexible. Lower oxygen permeability than rubber, aging resistant and non-oxidizing.

Quick and easy to plug onto fittings.

## Additional Information

---

No.	2-4461
Manufacturer (Brand)	neoLab
EAN	4058072003524
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polyvinyl chloride (PVC. Tygon)
Form	smooth straight
Length	15 m
DM outside	12.7 mm
DM inside	9.5 mm
TBST MAX	74 °C
TBST MIN	-50 °C
Wall thickness hose	1.6 mm
Application hose	for general laboratory applications
Hardness	55 shore A
for medium	Liquids Aqueous solutions Air/Gases/Vapors

