



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Tygon tubing S-hardness 55, 3.2 x 6.4 mm, 30 m/roll

Special Price
€39.00 was
€82.20
plus VAT &
Shipping

Product Images



Description

Ideal for general laboratory use.

- high chemical resistance
- clear, soft and flexible
- Shore hardness 55
- lower oxygen permeability than rubber
- resistant to aging and non-oxidizing
- thermally resistant from -50°C to +74°C

Additional Information

No.	2-4462
Manufacturer (Brand)	neoLab
EAN	4058072168643
Transport temperature	Room temperature
autoclavable at 121°C	No
Color	transparent (clear)
Material	Polyvinyl chloride (PVC. Tygon)
sterile	No
Form	smooth straight
Length	30 m
DM outside	6.4 mm
DM inside	3.2 mm
TBST MAX	74 °C
TBST MIN	-50 °C
Application hose	for general laboratory applications
Hardness	55 shore A
for medium	Liquids Aqueous solutions Air/Gases/Vapors

