



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Tygon silicone tubing USP VI/FDA 2.4 x 4.0 mm, 15 m/roll

€111.40
plus VAT &
Shipping

Product Images



Description

The right choice for laboratory, biotechnology and pharmaceutical applications!

Highly elastic with good restoring force for consistent fluid transfer, even with intermittent pump application. Ozone and UV resistant.

Particularly smooth inner surface.

PU: 15 m/roll

Additional Information

No.	2-4442
Manufacturer (Brand)	neoLab
EAN	4058072003340
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polyvinyl chloride (PVC. Tygon) Silicone
Form	smooth straight
Length	15 m
Hoses properties	food safe
Application hose	for general laboratory applications for peristaltic pumps
DM outside	4 mm
DM inside	2.4 mm
Hardness	50 shore A
TBST MAX	204 °C
TBST MIN	-80 °C
Wall thickness hose	0.8 mm
for medium	Liquids Aqueous solutions Air/Gases/Vapors

