

Description

The right choice for laboratory, biotechnology and pharmaceutical applications!

Highly elastic with good restoring force for consistent fluid transfer, even with intermittent pump application. Ozone and UV resistant.

Particularly smooth inner surface.

PU: 15 m/roll

Additional Information

No.	2-4445
Manufacturer (Brand)	neoLab
EAN	4058072003371
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polyvinyl chloride (PVC. Tygon) Silicone
Form	smooth straight
Length	15 m
DM outside	8 mm
DM inside	4.8 mm
TBST MAX	204 °C
TBST MIN	-80 °C
Wall thickness hose	1.6 mm
Hoses properties	food safe
Application hose	for general laboratory applications for peristaltic pumps
Hardness	50 shore A
for medium	Liquids Aqueous solutions Air/Gases/Vapors

