



neolab Migge GmbH
 Rischerstr. 7-9
 69123 Heidelberg
 Deutschland
 +49 (0)6221 /
 8442-44
<https://www.neolab.de>
 e

Umsatzsteuer-
 Identifikationsnummer
 :
 DE 143 450 657



neolab® Tygon silicone tubing USP VI/FDA 8.0 x 12.7 mm, 15 m/roll

€411.10
 plus VAT &
 Shipping

Product Images



Description

The right choice for laboratory, biotechnology and pharmaceutical applications!

Highly elastic with good restoring force for consistent fluid transfer, even with intermittent pump application. Ozone and UV resistant.

Particularly smooth inner surface.

PU: 15 m/roll

Additional Information

No.	2-4448
Manufacturer (Brand)	neoLab
EAN	4058072003401
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Polyvinyl chloride (PVC. Tygon) Silicone
Form	smooth straight
Length	15 m
DM outside	12.7 mm
DM inside	8 mm
TBST MAX	204 °C
TBST MIN	-80 °C
Wall thickness hose	2.4 mm
Hoses properties	food safe
Application hose	for general laboratory applications for peristaltic pumps
Hardness	50 shore A
for medium	Liquids Aqueous solutions Air/Gases/Vapors

