




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex Jet Pen Extension Tubing, 60 cm for Calibrex™ 525/530 25 ml

€98.90
plus VAT &
Shipping

Product Images



Description

Extension tubing for Calibrex™ dispenser

The Jet-Pen ejection cannula made of PTFE and the expandable spiral tubing provide maximum comfort when dispensing into vessels, up to a distance of 60 cm.

Easy assembly; completely autoclavable.

Additional Information

No.	SR-0248
Manufacturer (Brand)	Socorex
old neoLab article no.	7-7625
EAN	4058072367664
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (milky)
Material	Perfluoroethylene propylene (FEP)
Length	600 mm
DM inside	4 mm
Hoses properties	with special chemical resistance
Application hose	for Socorex Dispenser
for medium	Aqueous solutions

